



Brewer
more possibilities

BASIC EXAM TABLE

The Brewer Basic Exam Table

is clearly the:

- most reliable
 - most patient-friendly
 - most user-friendly
- entry exam table in the industry.

See why Brewer means... more possibilities.



Durability

- 3-year warranty

Safety

- Sturdy, retractable patient step with mat
- 325 lb. weight capacity*

Function

- Simple pneumatic back for easy patient positioning
- Easy to clean seamless upholstery
- Largest storage capacity available
- Front and side drawers
- Leg extension
- Plastic debris tray
- Pelvic tilt capability (optional)

Patient Comfort

- Sturdy step surface with mat
- Multi-density, bolstered support upholstery
- Adjustable ergonomic stirrups

Weight Capacity

- 325 lbs.*

Table Dimensions

- Height: 32"
- Length: 57.5" extends to 71"
- Width: 26.5"
- Leg extension length: 15"

Aesthetics

- 18 standard upholstery color options
- Custom colors available

Paper Roll Holder

- Accepts 21" x 3.5" and 18" x 3.5"

Storage

- Front drawers: 19.7" x 19.2" x 4"
- Side drawers: 22.7" x 20.2" x 7"
- Side drawers available on right or left side

Certifications

- UL listed, CUL listed

Electrical Models

- 120 VAC, 60 Hz, 2 AMPS

Options

- Hospital-grade duplex electrical outlet available

#4000
#4001 with drawer heater and pelvic tilt

accessories
#99505 stainless steel treatment pan
(12.75" x 10.38" x 2.5")
#100407 front drawer dividers

*under normal use

All Brewer Basic Exam Tables come with a 3-year warranty.

The Brewer Company reserves the right to improve or modify the design of its products consistent with new methods and techniques.

Medical Device Depot
Inc.

www.medicaldevicedepot.com
1-877-646-3300

Brewer
more possibilities