

FEATHERLIGHT

SERIES Delivers unequal flexibility with one-hand, leading counter-balanced movement.



Dichroic film-coated reflectors filter out 80% of beam heat

Three color-corrected, low voltage, 50-Watt halogen bulbs w/ heat sink technology that keeps the housing cool

"No-drift," counter-balanced arm design; 540° rotation at both the light head and yoke
Arm length: 36.6"

- Ideal for dermatology, plastic surgery, and minor procedures in the clinical environment.
- Bright, crisp light with no heat
- Movement provides effortless and accurate light placement without drifting or repositioning
- 5000 Footcandles @ 24"
- Color temperature - 4000° Kelvin
- Fused for your protection
- Bulb life: 4000 hrs; light pattern: 10" @ 24"
- Floor/Wall/Ceiling models
- Removable and sterilizable handle
- UL certified & CAN/CSA-C22.2 No.601.1 certified
- Ceiling models maximum length: 8': 25"; 9': 37"; 10': 49"

Diameter: 25"
Height: 5"
Compact five-leg caster base; power cord exits the base at floor level

Medical Device Depot Inc.
www.medicaldevicedepot.com
1-877-646-3300

accessories

#016150 50-watt replacement bulb

