

# DIGITAL HOLTER BY MIDMARK



**IQ**  
holter



The Midmark IQholter™ System combines a compact digital recorder with three powerful software applications. Our three tiers of Holter solutions are specifically designed to meet different practice needs, from family physicians to cardiologists.

# A POWERFUL HOLTER SYSTEM DESIGNED FOR YOUR PRACTICE.

## **IQholter™ Recorder:**

- Small and lightweight, the recorder weighs only 2.8 oz.
- Well-lit display and intuitive keypad make clinician navigation and menu control easy
- On-screen leads and channel identification for speed and simplicity during patient hookup
- Automatic startup self-test with warning alarm ensures reliable recording
- 5-lead hookup for 3-channel recording; 7-lead option available
- 24 hour standard recording (48 and 72 hour recording capabilities available)
- Removable Secure Digital (SD) card allows for immediate patient turnaround

## **IQholter™:**

- Color-coded arrhythmia events to identify different ECG morphologies
- Comprehensive reports with auto narrative summary to effectively reveal important test data and statistics
- Ability to edit and review test results and full disclosure locally or on the network without printing

## **IQholter™ EX:**

- Provides the ability to preview and adjust criteria settings to match patient's recording before the analysis
- Identifies Atrial Fibrillation (AF) events based on parameters set by the clinician

- Smart Page Scan to highlight specified events and reveal entire episodes at a glance
- Arrhythmia Templates with Intelligent Sorting to facilitate quick reviews and edits
- QT and QTc analysis for identifying prolonged intervals

## **IQholter™ EP:**

- Offers all the standard features of IQholter™ and IQholter™ EX
- Includes pacemaker detection to identify single and dual-chamber paced beats for non-captures and sensing failures
- Automatic pacemaker detection eliminates setting changes between patients



The IQholter™ Recorder can be worn as comfortably as an MP3 player or cell phone.

# IQholter™ Technology by Midmark

All the features you need, in a solution tailored just for you.

Designed with a well lit display and intuitive keypad for ease of use.



Easy to use SD card allows for immediate patient turnaround.



Comprehensive report with auto narrative summary.



Fully automated analysis and reporting process.



Make your practice more efficient with an IQ Digital Diagnostic System  
**Workflow efficiency and improvements in patient care are just a few of the benefits that can result from implementing an IQ Digital Diagnostic System. Call us today and we'll put together a customized suite of solutions to help grow as your practice grows.**

Diagnostic Products									
	ECG	Holter	Stress	Spirometer	Vital Signs	EHR/EMR Connectivity	Single Platform (IQmanager™)	Standalone	Thin - Client (IQpath™)
IQecg™	■					■	■		■
IQflex™	■							■	
IQclassic™	■							■	
IQholter™		■				■	■		■
IQholter™ EX		■				■	■		■
IQholter™ EP		■				■	■		■
IQstress™			■			■	■		
IQspiro™				■		■	■		■
IQvitals™					■	■	■	■	■
IQvitals™ PC					■	■	■		■

### Minimum System Requirements

- Windows® Vista Home, Windows® XP Pro, Windows® 7 Home Premium, 32-bit and 64-bit
- Intel® Core™ 2 Duo Processor E4300 or 64-bit processor or faster
- 2 GB of RAM
- 2 GB available hard disk space
- USB or Serial Port
- Mouse, keyboard and high resolution display monitor

Please see the IQholter™ Operation Manual or visit [midmark.com](http://midmark.com) for more detailed system requirements

### Ordering Information

4-000-0110	IQholter™ Digital Holter with Recorder	1-370-0020	IQholter™ Upgrade Kit (24, 48, 72 Hours)
4-000-0113	IQholter™ EX Digital Holter with Recorder	2-100-0080	IQholter™ Prep Kit, 5-lead
4-000-0116	IQholter™ EP Digital Holter with Recorder	2-100-0090	IQholter™ Prep Kit, 7-lead
4-000-0119	IQholter™ Digital Holter Recorder (5-lead standard)	2-100-0060	IQholter™ Electrodes (box of 50)
		2-100-0050	IQholter™ Recorder Pouches, Disposable (pack of 25)
4-100-1200	IQmanager™ Software - Patient data management of IQecg™, IQspiro™, IQholter™, IQstress™, IQvitals™, and IQvitals™ PC	2-100-0051	IQholter™ Recorder Pouch with Belt, Reusable
1-370-0010	IQholter™ Secure Digital Memory Card	90-30-0205	IQholter™ Lead Wire, 5-lead
		90-30-0206	IQholter™ Lead Wire, 7-lead



## RETHINK THE CLINICAL SPACE

At the core of Midmark is a new way to think about caregiver/patient interaction and the delivery of efficient patient care. We are dedicated to working together with you to redefine the future of the clinical space.

[midmarkclinicalsolutions.com](http://midmarkclinicalsolutions.com)

For information, pricing or to place an order, contact us at: