

**NEW!
KEYPAD FEATURE!**

LifeDop™ 300 ABI

by SummitDoppler

**Economical Hand-Held Doppler System with
FASTER ABI Calculation & Waveform Printout**



Perform the ABI exam for P.A.D. diagnosis
FASTER & EASIER with the LifeDop 300 ABI!

- Keypad for easily entering pressure data & calculating the indices
- Stores pressure data & two ankle waveforms prior to printing
- Proven technology: bi-directional probe and high speed printer
- Improved flow direction indicator on display

**QUICK
ROI!
\$121 per
exam***

To order the LIFEDop 300 ABI use:

Model # L300AC: ABI System with 4 Cuffs & Aneroid

Medical Device Depot
Inc.

www.medicaldepositedot.com
1-877-646-3300

LifeDop™ 300 ABI

Economical Hand-Held Doppler System with FASTER ABI Calculation & Waveform Printout



ABI 4 Cuff Package includes:

- Doppler main unit
- 8 MHz bi-directional probe
- Portable printer with cable
- Aneroid
- (2) 10 cm cuffs
- (2) 12 cm cuffs
- Adhesive-backed label paper
- Report forms
- ABI training video
- ABI chart
- Ultrasound gel (60 g)
- Carrying case

SummitDoppler

**QUICK
ROI!**
\$121 per
exam*

To order the LifeDop 300 ABI use:

Model # L300AC: ABI System with 4 Cuffs & Aneroid

Main Unit and Probe

- Bi-directional 8 MHz high sensitivity probe (optional 5 MHz probe)
- Separate waveforms indicate flow toward and away from the probe
- Auto-scaling and no calibration required
- 200-4,000 Hz bandwidth covers peripheral applications
- User replaceable probe cable
- Low battery indicator
- Long battery life (3 AA alkaline batteries)
- Dimensions: 147 x 70 x 35 mm (h w l)
- Weight: 320 grams

ABI Calculation and Waveform Storage

- Calculates ABI for both sides using the higher of the two arm pressures
- Accepts pressures from 30 to 299 mmHg
- Stores Doppler waveform from each ankle for printing

Compact Portable Printer

- Label paper for fast report generation
- High speed
- Rechargeable (Li-Ion battery)
- Size: 136 x 84 x 43 mm (h w l)
- Weight: 280 grams



**Perform the ABI exam for P.A.D. diagnosis
FASTER & EASIER with the LifeDop 300 ABI!**

Medical Device Depot Inc.

www.medicaldevicedepot.com

1-877-646-3300