

LL100™ Cryosurgical System

gives physicians the fastest freeze and defrost

3525H Ellicott Mills Drive, Suite 108
 Ellicott City, MD 21043
 877-646-3300 phone
 410-505-4567 fax



Wallach's exclusive two-trigger design lets the physician precisely control both freeze and defrost. With twice the freeze power of other systems, the LL100 performs both functions more effectively in less time eliminating treatment failures. Freeze and defrost only occur when triggered. Light weight, well-balanced and maneuverable – provides full view of the target tissue. The gas hose allows for unimpeded surgical movement. Cryotips can be changed at any time without venting. A large selection of highly conductive reusable, sterilizable cryotips and disposable plastic shields are available.

REF 900001	LL100 (N₂O)
REF 900002	LL100 (CO₂)
REF 901060	N₂O Cylinder (unfilled)
REF 901061	CO₂ Cylinder (unfilled)
REF 901011	Cart for 20 lb. Cylinder

5 year warranty*

*see Warranty Policy on page 32

WA1000B Cryosurgical Console System

*High performance features and professional appearance
 make the WA1000B your best choice*

The WA1000B console offers the ideal combination of easy-to-use controls and built-in safety features. It uses the powerful LL100 two-trigger freezer. The Temperature Indicator with integrated timer gives physicians precise documentation during cryosurgery. In addition, a removable stainless steel well holds the freezer and may be used as a disinfecting receptacle. The cart can carry a single 20 lb. cylinder or twin "E" cylinders.

REF 900506-1	WA1000B with Blank Console, Base and LL100 Multi-Tip Freezer
REF 900506-2	WA1000B with Digital Temperature Indicator, Timer, Base and LL100 Multi-Tip Freezer with Thermocouple

Cryosurgical Tips

highly conductive, thermal, reusable, sterilizable cryotips with disposable plastic shields

Medical Device Depot^{Inc.}

- Utilized by LL100 and WA1000B Cryosurgical Systems
- Ultimate in safety and convenience
- Each tip is supplied with 5 single-use, disposable plastic shields
- Remove disposable shield from cryotip and discard after each use (shields insulate the tip and assure proper attachment)
- Cryotip is fully sterilizable to prevent cross-contamination (autoclave or ETO)

www.medicaldevice depot.com

3525H Ellicott Mills Drive, Suite 108
Ellicott City, MD 21043
877-646-3300 phone
410-505-4567 fax



2.5mm HPV Tip (T-0200)
REF 900300AA



Exocervical (T-1905)
REF 900209AA



Procto Curved (T-0832)
REF 900217AA



5.0mm HPV Tip (T-0500)
REF 900301AA



Exo-Endocervical (T-1910)
REF 900211AA



Anorectal (T-0337)
REF 900218 (non-sterilizable)



7.0mm HPV Tip (T-0700)
REF 900302AA



Long Exo-Endocervical (T-1920)
REF 900212AA



Anorectal (T-0337)
REF 900220 (non-sterilizable)



9.0mm HPV Tip (T-0900)
REF 900303AA



Exo-Endocervical (T-1920)
REF 900210AA



Bevel (T-0524)
REF 900202AA



11.7mm HPV Tip (T-1200)
REF 900304AA



Exocervical (T-2507)
REF 900214AA



Bevel (T-0823)
REF 900207AA



13.5mm HPV Tip (T-1300)
REF 900305AA



Exo-Endocervical (T-2514)
REF 900215AA



Cone (T-0219)
REF 900201AA



15.0mm HPV Tip (T-1500)
REF 900306AA



Condyloma (T-0519)
REF 900203AA



Microderm (T0219)
REF 900200AA



Flat (T-1900)
REF 900208AA



General Purpose (T-0819)
REF 900205AA



Cone (T0826)
REF 900206AA



Flat (T-2500)
REF 900213AA



Procto (T-0832)
REF 900204AA

REF 900221AA (Bag of 100)

Disposable Plastic Shields

\$79.00

