



Propaq LT | Highly Deployable Continuous Monitor

Introducing the Innovative New  
**Welch Allyn Propaq® LT**

It won't just change how you monitor patients.  
It'll change where you do it.



**Medical Device Depot**  
Inc.  
[www.medicaldevicedepot.com](http://www.medicaldevicedepot.com)  
1-877-646-3300

**Welch Allyn®**

Advancing Frontline Care™

# The New Welch Allyn Propaq® LT

Monitor more patients in more places than ever before.



Reliability, durability and affordability combine to give you perhaps the most versatile monitor on today's market.

The new Propaq LT is the perfect monitor for today's fast-paced frontline hospital environment—letting you extend continuous monitoring to areas in the hospital where technical complexity or cost made it impossible in the past. Because it's so easy to use and configure, the new Propaq LT is the perfect choice for nurses on general care floors. And with legendary Propaq quality, you know it's rugged and reliable enough to use in the ED and for interdepartmental transport.

With wireless connectivity and spot-check mode options, the new Propaq LT is like getting three monitors—ambulatory, bedside and transport—for the price of one.



The one monitor you can deploy when—and where—you need it most.



Equipped with the optional Propaq LT large color display and charging cradle, the Propaq LT is ideal for use as a bedside or procedure monitor. And with a revolutionary new design that's light enough to fit on an IV pole, rugged enough for air transport, and versatile enough for pediatric applications, the Propaq LT gives you the power to monitor more patients in more places than ever before.

SHOWN ACTUAL SIZE:  
The Welch Allyn Propaq LT  
weighs about two pounds.



Tab/Wave	80	127/71	12	97
Time	HR/min	NIBP mmHg (86)	Resp/min	SpO2
12:30	80	n/a	12	97
12:29	80	n/a	12	97
12:28	80	n/a	12	97
12:27	80	n/a	12	97
12:26	80	127/71 (86)	12	97
12:25	80	n/a	12	97
12:24	80	122/73 (86)	12	97

## Key Features:

- Shock resistance of up to 75 g—enough to withstand a 6-foot drop
- Heart rate, 3- or 5-lead ECG, SpO<sub>2</sub>, NIBP, respiration
- Neonatal, pediatric and adult patient modes
- Continuous monitoring and spot-check modes
- 24-hour tabular and graphical trending
- Storage of up to 20 full-disclosure data snapshots
- Wireless connectivity option to Welch Allyn Acuity® Central Station via FlexNet™
- Designed for use in large- and small-body fixed- and rotary-wing aircraft



Quickly toggle between three display modes that are right for your care area (five to choose from). Shown are two-waveform, tabular trend with waveform, and large numeric (pictured above) displays.