

BEAM

Omni Lifter Stand

The unique OmniLifter Irrigation stand is the answer for your fluid-intensive IV irrigation procedures. This stand uses a foot pump to easily raise and lower solutions ensuring continuous fluid delivery during surgery. The ergonomic foot pump system greatly reduces the threat of back strain or injury related to improper lifting.

Features include:

- Sturdy, stable base with two locking casters provides “tip-free” usage
- Foot pump continuously lifts and lowers up to 6-3,000 cc or 4-5,000 cc bags
- Pole decent is controlled by an adjustable valve allowing smooth and safe operation.
- Unit constructed of heavy-duty aluminum with stain resistant finish for longevity and sanitation
- Pole locks in position with a clamp for safe operation

Free demonstration and 30-day trial available

	Heavy Weight™ Infusion Stand #741300	Jr. Heavy Weight™ Infusion Stand #741301	Arthroscopy Stand #741302	Lite Weight Stand #741303	Omni Lifter Stand #741313
Dimension	Max Height 9'	Max Height 7'	Max Height 8½'	Max Height 7'	Max Height 9'
Price					

NEW!



Medical Device Depot Inc.
www.medicaldecedepot.com
 Toll Free 1-877-646-3300



#741313