

ECHO IMAGE SERIES

Medical Device Depot Inc.

www.medicaldevicedepot.com
Toll Free 1-877-646-3300

FEATURES

- 2" high-density urethane foam padding
- The top and frame are upholstered together to prevent dirt from collecting in seams
- Stationary models feature easy to clean laminate surfaces available in several colors
- All metal slide bar supports window drop
- 3 bariatrics models on page 18
- Power units feature:
 - All welded steel frame
 - Single drive motor system
 - Hospital-grade 110 volt plug
 - Fully enclosed lift system
 - Hands-free foot control

	Length	Height	Width
4615	76"	31"	27"

Echo Image Table

- All laminate base
- Multi-position adjustable headrest
- Full triple leg support
- 10" wide x 12" deep drop leaf window
- 400 lbs. capacity



Double Window Drop Supports



Top and Contoured Frame Upholstered Together



	Length	Height Range	Width
4620	76"	20"-30"	33"

Power Echo Image Table

- All welded steel base and frame
- Dependable, safe, single motor, enclosed lift mechanism
- Multi-position headrest
- Easy wheelchair access
- Hospital-grade plug
- Hands-free foot control
- 8" wide x 11" deep drop leaf window
- 110 volt/60 hz./2.8 amps
- 400 lbs. capacity



	Length	Height Range	Width
4630	76"	20"-30"	33"

Power Echo Image Table


- All welded steel base and frame
- Dependable, safe, single motor, enclosed lift mechanism
- Multi-position headrest
- Easy wheelchair access
- Hospital-grade plug
- Hands-free foot control
- Built-in flip-down casters for easy movement when not in use
- 8" wide x 11" deep drop leaf window
- 110 volt/60 hz./2.8 amps
- 400 lbs. capacity



	Length	Height Range	Width
4690V	71"	24"-34"	32"

Vascular Table

- 2.5" thick, high density, foam padding
- All welded steel base and frame
- Fold-down echo window
- Electric adjustable height range
- Electric tilting table top (head and foot end tilt 20°)
- Protractor for table top angle measurement
- Manual lift back with switch activated gas spring
- 8" wide x 11" deep fold-down section
- Hands-free foot controls
- 110 volt/60 hz./2.8 amps
- 400 lbs. capacity
- Approximate weight 225 lbs..
- Choice of upholstery colors with grey base



Tilting Positions