

ETA TREATMENT CABINETS

Clinton ETA (Easy-To-Assemble) Treatment Cabinets are designed with the needs of the professional in mind. Six models feature modules that ship easy and go together fast using only a Phillips screwdriver. Mix or match up to 12 colors of laminate for a traditional or unique look.

SHIPPING ADVANTAGES

- Save money on shipping. Clinton ETA cabinets ship via UPS® or FedEx® Ground
- Single person delivery and assembly
- No inside delivery charges

NOTE:

Large orders that are less expensive to be shipped by common carrier will be shipped assembled.

FEATURES


- Pre-assembled modules
- All assembly holes pre-drilled and pre-threaded
- Easy-clean all laminate construction
- Tough chip-resistant PVC edges
- Modules are assembled with steel fasteners, not wood dowels
- 3-way, adjustable 110° door hinges
- Euro-style drawers with metal sides that hold up to 75 lbs.
- Recessed toe kick or dual wheel caster bases
- Height-adjustable drawer fronts
- Mix color choices for top sides and fronts
- Light color interior laminates
- Any damaged components are field changeable
- Assemble with only a Phillips screwdriver

ORDERING INFORMATION:

When ordering cabinets you need to select top, front and side colors. If only one color is selected, that color will be used for the entire cabinet.

Medical Device Depot
Inc.

www.medicaldevicedepot.com
1-877-646-3300



Clinton uses pre-drilled holes and pre-threaded fasteners that allow Clinton ETA cabinets to be assembled faster and easier.



	Length	Height	Depth
8821	24"	35"	18"

Treatment Cabinet with 2 Doors & 1 Drawer

- 1 adjustable shelf
- Euro-style drawer with metal sides that hold up to 75 lbs.
- All laminate construction

	Length	Height	Depth
8822	24"	35"	18"

Treatment Cabinet with 2 Doors & 2 Drawers

- 1 adjustable shelf
- Euro-style drawers with metal sides that hold up to 75 lbs.
- All laminate construction

	Length	Height	Depth
8805	24"	35"	18"

Treatment Cabinet with 5 Drawers

- Smooth glide tracks for easy open and close
- Euro-style drawers with metal sides that hold up to 75 lbs.
- All laminate construction



	Length	Height	Depth
8922	24"	34"	18"

Mobile Treatment Cabinet with 2 Doors and 2 Drawers

- Mounted on dual wheel swivel casters
- Euro-style drawers with metal sides that hold up to 75 lbs.
- All laminate construction

	Length	Height	Depth
8921	24"	34"	18"

Mobile Treatment Cabinet with 2 Doors and 1 Drawer

- Mounted on dual wheel swivel casters
- Euro-style drawer with metal sides that hold up to 75 lbs.
- All laminate construction

	Length	Height	Depth
8950	24"	34"	18"

Mobile Treatment Cabinet with 5 Drawers

- Mounted on dual wheel swivel casters
- Euro-style drawers with metal sides that hold up to 75 lbs.
- All laminate construction