

FUN SERIES TABLES

TREATMENT TABLE FEATURES

- All graphics are permanently embedded in the high pressure laminate and will not peel off
- Paper dispenser included
- Optional locking doors available
- Light color interiors
- All colors are as shown

SCALE TABLE FEATURES

- Comfortable work height of 36"
- Top is rimmed on three sides for safety
- Top's front and back rim edges are beveled and the corners rounded for safety
- Easy-clean, replaceable center pad
- Built-in infantometer with folding ends for measuring height
- Easy access to weight, measurement and storage
- Built-in concealed paper dispenser and cutter
- Recessed front toe kick
- Adjustable built-in levelers provide a stable surface
- Easy-clean all laminate base
- Optional storage area locks available
- Ample 2-door storage
- 3 innovative tables to choose from

SCALE FEATURES


- Dependable, accurate Tanita scale
- Easy-read digital display
- Auto calibration on start-up
- Auto shut-off
- Electric or battery operation (6 AA batteries required)
- Weights in increments of lbs., ozs. or kilos
- Maximum weight 40 lbs.



Top Rimmed on 3 Sides



Beveled Sides



Infantometer



Digital Readout

Medical Device Depot Inc.

www.medicaldevicedepot.com
1-877-646-3300

FUN SERIES SCALE TABLES



	Length	Height	Width
7822 top	43"	–	24"
base	47"	36"	21"



Scale Table/Zoo Bus with Jungle Friends

- Features 1 end door and 1 side door with 1 adjustable shelf



	Length	Height	Width
7833 top	43"	–	24"
base	47"	36"	21"



Scale Table/Engine K-9 with Dalmatian Firefighters

- Features 1 end door and 1 side door with 1 adjustable shelf



	Length	Height	Width
7844 top	43"	–	24"
base	47"	36"	21"



Scale Table/Outback Buggy with Aussie Animal Pals

- Features 1 end door and 1 side door with 1 adjustable shelf


FUN SERIES TREATMENT TABLES



	Length	Height	Width
7020	72"	34"	27"

Zoo Bus with Jungle Friends

- Features 4-door storage with 3 adjustable shelves



	Length	Height	Width
7030	72"	34"	27"

Engine K-9 with Dalmatian Firefighters

- Features 4-door storage with 3 adjustable shelves



	Length	Height	Width
7040 top	43"	-	27"
base	47"	35"	21"

Outback Buggy with Aussie Animal Pals

- Features 2-door storage, removable center pad and top rimmed with stainless steel tubing on 3 sides



Medical Device Depot Inc.


www.medicaldevicedepot.com
1-877-646-3300



	Length	Height	Width
7050	72"	34"	27"

Cool Camper

- Features 4-door storage with 3 adjustable shelves



	Length	Height	Width
7060	68"	34"	27"

Medicine Show Wagon

- Features 3-door storage with 3 adjustable shelves