

# Hill Adjustable HA90

*The Classic  
Variable-Height  
Table Design*

**Medical Device Depot**  
Inc.

www.medicaldevicedepot.com  
1-877-646-3300



Pictured with the optional face cut-out, the HA90 may also be ordered with a gas spring-assisted lift-back that raises with ease to any position up to 80°. An electric lift-back with foot pedal is also available.

During manipulation, therapy, or mobilization, the HA90 makes your job easier by elevating to the perfect work height. You will appreciate the reduction in

bending and back stress and your patients will appreciate easy access onto the table provided by the variable height. Gliding quietly and smoothly, the HA90 elevates from 20" to 30" with a touch of the foot pedal. All mechanisms are completely enclosed to ensure safety for you and your patients.

- Standard 24" wide top
- Standard 1 1/2" foam
- Optional widths of 27" and 30"
- Optional fold-under stirrups
- Optional extra foot pedal,
- Optional face cut-out, with or without insert
- Optional Side Rails

## HA90 Options



Rocker foot pedals can be substituted for the corded pedal and are mounted on both sides of the table.



Casters add mobility to your Hill table and are engaged with the foot levers from either side of the table.



A storage shelf can be added to the head or foot end of the table for convenient access to often used items.



A crescent headpiece is available for massage therapy and includes soft foam and adjustable lock positioning.