

# For Aesthetics

## HA90ALB

The 90ALB (Aesthetic Lift Back) is an extremely affordable chair for skin care with 2" of soft foam over 1 1/2" of medium-density for maximum patient comfort. Many standard and optional features make the 90ALB ideal for aesthetic work.



- Shown with  
Bordeaux upholstery  
and Beige base.


### Easy up. Easy down.

The 90ALB features electric elevation from 23" to 36" or 21" to 31". The standard model has an assisted, gas spring back for easy positioning, an adjustable headrest, fold-down removable armrests, contour back and durable, upgraded designer upholstery.

The HA90ALB can fully recline and is also available with an optional electric lift back and a variety of head cushion options.



# Specifications



Powered by an electric actuator 60 Hz, 110 volts

**Amperage:**

One Motor Table - 2.4, Two Motor Table - 5.1  
 Three Motor Table - 7.5, Four Motor Table - 10.2

**Lifting Capacity:**

400 lbs. (600 lbs. capacity upon request)

**Height Range:**

90ALB - 23" to 36" or 21" to 31"  
 90D - 22 to 35. With Trendelenburg - 24" to 37"  
 or optional 22" to 32"

**Table Base:**

Beige or Gray telescoping base.  
 All metal parts are powder coated.

**Shipping Weight:**

90ALB - 225 lbs., 90D - 300 lbs.

**Warranty:** One year

## Color Choices



Upholster your table in durable Boltaflex® with Prefix™, scientifically formulated to give superior protection from staining and abrasions (top row). Or, for more of an upscale look and feel, choose puncture-proof Finesse™ material (bottom row). Custom colors are available at additional cost. (Shades may vary from actual material. Colors may change without notice.)