

INSTRUMENT STANDS

Medical Device Depot Inc.
www.medicaldevicedepot.com
 1-877-646-3300



	Base Diameter	Height Range
M-20	23"	31"–50"

Mobile Instrument Stand

- 5-leg design aluminum base
- Removable stainless steel tray
- Dual wheel soft roll casters
- Chrome plated pole and frame

	Base Diameter	Height Range
M-25	23"	31"–50"

Mobile Instrument Stand

- 5-leg design nylon base
- Removable stainless steel tray
- Dual wheel casters
- Chrome plated pole and frame



	Base	Height Range
M-22	21" x 20"	34"–54"

Double Post Mayo Stand

- Removable stainless steel tray
- 2 ball bearing, rubber casters
- 1" chrome plated poles and base
- Easy position base

	Base	Height Range
M-21	18" x 22"	31"–50"

Single Post Mayo Stand

- Removable stainless steel tray
- 2 ball bearing, rubber casters
- Chrome plated pole and base
- U-shaped base



	Length	Width	Height
M-24	19"	12 ³ / ₄ "	3/4"

Stainless Steel Replacement Trays

- Rimmed tray with satin finish



	Base Diameter	Height Range
M-26	22"	30"–49"

Mobile Instrument Stand

- 4-leg heavy-duty design
- Removable stainless steel tray
- Rubber, ball bearing casters
- Chrome plated pole and frame
- Powder coated base with bumper ends