

# INSTRUMENT STANDS

**Medical Device Depot** Inc.

www.medicaldevicedepot.com  
Toll Free 1-877-646-3300



|             | Base Diameter | Height Range |
|-------------|---------------|--------------|
| <b>M-20</b> | 23"           | 31"-50"      |

**Mobile Instrument Stand**

- 5-leg design aluminum base
- Removable stainless steel tray
- Dual wheel soft roll casters
- Chrome plated pole and frame

|             | Base Diameter | Height Range |
|-------------|---------------|--------------|
| <b>M-25</b> | 23"           | 31"-50"      |

**Mobile Instrument Stand**

- 5-leg design nylon base
- Removable stainless steel tray
- Dual wheel casters
- Chrome plated pole and frame



|             | Length | Width                            | Height |
|-------------|--------|----------------------------------|--------|
| <b>M-24</b> | 19"    | 12 <sup>3</sup> / <sub>4</sub> " | 3/4"   |

**Stainless Steel Replacement Trays**

- Rimmed tray with satin finish



|             | Base      | Height Range |
|-------------|-----------|--------------|
| <b>M-22</b> | 21" x 20" | 34"-54"      |

**Double Post Mayo Stand**

- Removable stainless steel tray
- 2 ball bearing, rubber casters
- 1" chrome plated poles and base
- Easy position base

|             | Base      | Height Range |
|-------------|-----------|--------------|
| <b>M-21</b> | 18" x 22" | 31"-50"      |

**Single Post Mayo Stand**

- Removable stainless steel tray
- 2 ball bearing, rubber casters
- Chrome plated pole and base
- U-shaped base



|             | Base Diameter | Height Range |
|-------------|---------------|--------------|
| <b>M-26</b> | 22"           | 30"-49"      |

**Mobile Instrument Stand**

- 4-leg heavy-duty design
- Removable stainless steel tray
- Rubber, ball bearing casters
- Chrome plated pole and frame
- Powder coated base with bumper ends