


# CLASSIC ECG BY MIDMARK



IQ  
classic





An ECG is often the first line of defense and a primary weapon against heart disease. Understanding the different options available with an ECG and how it fits into your office workflow can make for a daunting task. To make it easier, Midmark now offers a traditional ECG with best in class features that seamlessly fits into your daily workflow.

# TRUSTED TECHNOLOGY TO PROTECT AGAINST HEART DISEASE

## Simplicity

- Portable with an integrated printer, the IQclassic™ can be used in any exam room
- Perform ECG tests where you need to with internal battery operation
- Simple, one-touch buttons for quick and easy data entry, ECG acquisition, rhythm printing and transmission
- Backlit LCD display provides a real-time preview of 3-, 8- or 12-lead ECG
- Choose to automatically save up to 60 ECG tests with standard internal storage for retrieval and printing later

## Flexibility

- Choose from Multi Standard or Cabrera print options including: 3+1, 6, 12, 3+3 and 6+6 channels
- Choose from rhythm print options 3, 6 and 12 channels
- FAA modem option adds the ability to address FAA physicals for aviators
- Automatic transition between AC and battery operation allows clinicians to focus on patient care

## Reliability

- Adult and pediatric, 12-lead resting ECG interpretation algorithm by VERITAS™ is globally recognized as best in class in interpretation
- 10,000 samples per second per channel for accurate pacemaker spike detection
- 1,000 samples per second/ per channel for accurate ECG signal detection and analysis
- Provides reliable and accurate analysis with interpretation and reasoning statements for physicians



Portable with an integrated printer, the IQclassic™ can be used in any exam room

# IQclassic™ - An ECG that fits your needs

IQclassic™ is a 12-lead resting ECG, featuring the globally recognized, best in class VERITAS™ algorithm.

It is easy to use, flexible to configure and complements your workflow.



Standard full-size paper (8.5" x 11") output includes multiple standard or Cabrera print options.



Conduct FAA physicals for aviators with FAA modem capability.



Backlit LCD display provides a real-time preview of 3-, 8- or 12-lead ECG.



Full size keypad and dedicated one-touch buttons make it easy for staff to learn ECG acquisition and navigation.

Make your practice more efficient with an IQ Digital Diagnostic System  
**Workflow efficiency and improvements in patient care are just a few of the benefits that can result from implementing an IQ Digital Diagnostic System. Call us today and we'll put together a customized suite of solutions to help your practice evolve with changes in healthcare.**

Diagnostic Products				
	EHR/EMR Connectivity	Single Platform (IQmanager™)	Standalone	Thin - Client (IQpath™)
IQecg™	■	■		■
IQclassic™ ECG			■	
IQholter™	■	■		■
IQholter™ EX	■	■		■
IQholter™ EP	■	■		■
IQstress™	■	■		
IQspiro™	■	■		■
IQvitals™	■	■	■	■
IQvitals PC™	■	■		■

### Models

- 4-000-0070 IQclassic™, Serial RS232
- 4-000-0071 IQclassic™, Serial & Modem

### Accessories



Disposable EZ-Trode ECG Electrodes



Clear ECG Clips (4mm & Snap)



IQcart™



## RETHINK THE CLINICAL SPACE

At the core of Midmark is a new way to think about caregiver/patient interaction and the delivery of efficient patient care. We are dedicated to working together with you to redefine the future of the clinical space.

[midmarkclinicalolutions.com](http://midmarkclinicalolutions.com)

For information, pricing or to place an order, contact us at:

**Medical Device Depot**  
 Inc.  
[www.medicaldevicedepot.com](http://www.medicaldevicedepot.com)  
 1-877-646-3300