



	Base Diameter	Height Range
IV-30	22"	51 1/4"-93"

Value 2-Hook IV Pole

- 4-leg epoxy coated steel base
- Chrome plated steel pole
- Twist lock height adjustment
- Detachable ram's horn hooks with plastic tips
- Dual wheel casters
- Epoxy coated 4-leg base with bumper guard

	Base Diameter	Height Range
IV-31	22"	51 1/4"-94"

Heavy Base 2-Hook IV Pole

- 4-leg epoxy coated steel base
- Chrome plated steel pole
- Single knob height adjustment
- Detachable ram's horn hooks with plastic tips
- 2" rubber wheel ball bearing casters

	Base Diameter	Height Range
IV-314	22"	51 1/4"-94"

Heavy Base 4-Hook IV Pole

- 4-leg epoxy coated steel base
- Chrome plated steel pole
- Single knob height adjustment
- Detachable ram's horn hooks with plastic tips
- 2" rubber wheel ball bearing casters

	Base Diameter	Height Range
IV-35	24"	50"-90"

Five-Leg 2-Hook IV Pole

- 5-leg aluminum base
- Chrome plated steel pole
- Single knob height adjustment
- Detachable ram's horn hooks with plastic tips
- Soft roll, dual wheel casters

	Base Diameter	Height Range
IV-33	26"	56"-99"

Six-Leg 2-Hook IV Pole

- 6-leg epoxy coated base
- Large diameter chrome plated steel pole
- Single knob height adjustment
- Detachable ram's horn hooks with plastic tips
- 3" rubber wheel, ball bearing casters

	Base Diameter	Height Range
IV-37	24"	51"-91"

Five-Leg 2-Hook IV Pole

- 5-leg, easy-clean, nylon composite base
- Chrome plated steel pole
- Single knob height adjustment
- Detachable ram's horn hooks with plastic tips
- Dual wheel casters

IV STANDS



	Base Diameter	Height Range
IV-354	24"	50"–90"

Five-Leg 4-Hook IV Pole

- 5-leg aluminum base
- Chrome plated steel pole
- Single knob height adjustment
- Detachable ram's horn hooks with plastic tips
- Soft roll dual wheel casters

	Base Diameter	Height Range
IV-334	26"	56"–99"

Six-Leg 4-Hook IV Pole

- 6-leg epoxy coated base
- Large diameter chrome plated steel pole
- Single knob height adjustment
- Detachable ram's horn hooks with plastic tips
- 3" rubber wheel, ball bearing casters

	Base Diameter	Height Range
IV-374	24"	51"–91"

Five-Leg 4-Hook IV Pole

- 5-leg, easy-clean, nylon composite base
- Chrome plated steel pole
- Single knob height adjustment
- Detachable ram's horn hooks with plastic tips
- Dual wheel casters

IV STANDS COMPARISON CHART

Model No.	Hooks	Base	Base Diameter	Caster	Height Range	Load Capacity
IV-35	2	Aluminum	24"	5 Dual Wheel	50"–90"	45 lbs.
IV-354	4	Aluminum	24"	5 Dual Wheel	50"–90"	45 lbs.
IV-33	2	Epoxy/Steel	26"	6 Rubber	56"–99"	45 lbs.
IV-334	4	Epoxy/Steel	26"	6 Rubber	56"–99"	45 lbs.
IV-37	2	Nylon	24"	5 Dual Wheel	51"–91"	45 lbs.
IV-374	4	Nylon	24"	5 Dual Wheel	51"–91"	45 lbs.
IV-30	2	Epoxy/Steel	22"	4 Dual Wheel	51¼"–93"	—
IV-31	2	Epoxy/Steel	22"	4 Rubber	51¼"–94"	45 lbs.
IV-314	4	Epoxy/Steel	22"	4 Rubber	51¼"–94"	45 lbs.