

# MiniVID

The MiniVID eyepiece camera is available with an RCA output or USB connection. Plug it in to your TV, VCR, or computer. Then, attach it to any microscope for instant video micrography. The innovative MiniVID eyepiece camera slips into any standard microscope eyepiece tube, in place of the eyepiece, for crisp imaging in an instant.

## MINIVID USB

### Image Sensor

3.1M pixel CMOS chip

### Max Pixels

2048 X 1536 pixels

### Interface

USB 2.0

### Exposure

Automatic

### Saved Picture File Format

BMP, JPG

### System Requirement

Windows 98, XP, or Vista

### Software

Scopephoto Senior micro-image analysis software



## MiniVID

Model #	Description
MVC-U3MP-EMT3	MiniVID USB eyepiece video camera
MVC-RN77-EMT1	MiniVID RCA eyepiece video camera

## MINIVID RCA

### Image Sensor

1/3 inch color CMOS chip

### Output

6-foot cord with composite video (RCA)

### Sensitivity

5 lux

### Signal

NTSC or PAL available

### Power

9V DC < 1amp

Includes 110V AC adapter, 220V available



  
LW  
Scientific

**Medical Device Depot**  
Inc.

www.medicaldevicedepot.com  
1-877-646-3300