

Nicolet Pocket-Dop II

Obstetric and Vascular Versatility

Pocket-Dop II is a full-featured Doppler with interchangeable obstetric and vascular probes. Use the 2 or 3 MHz OB probes to monitor the fetal heartbeat. Use the 5 or 8 MHz vascular probes to screen patients at risk for stroke, peripheral arterial disease and other occlusive vascular conditions or to monitor difficult pulses and obtain systolic blood pressures quickly and easily.

Obstetric Examination

The Pocket-Dop II is a valuable aid in detecting pregnancy and monitoring the fetal heartbeat. The built-in speaker amplifies fetal sounds to reassure the expectant mother. A complete obstetric system should include a 3 MHz probe for early fetal heartbeat detection and a 2 MHz probe for use in later pregnancy or with larger patients.

A 2 MHz waterproof probe is ideal for underwater labor or delivery.

Vascular Applications

Use the Pocket-Dop II to compare simple Doppler blood pressures in the arm and ankle to help diagnose arterial disease in the lower extremities. Or quickly screen for lower extremity vascular conditions such as valvular incompetence.

Emergencies and Challenging Clinical Situations

The Pocket-Dop II can be used to accurately monitor pulses and systolic blood pressures, even in the noisiest environments; confirm blood flow to wound sites, locate vessels for catheterization, withdrawal, injection or IV therapy; clearly hear low systolic blood pressures with patients who are diabetic, elderly or in shock.

Catalog Numbers

Pocket-Dop II with:

2 MHz probe	P200
3 MHz probe	P300
5 MHz probe	P500
8 MHz probe	P800
2 & 3 MHz probes	P230
2 & 5 MHz probes	P250
2 & 8 MHz probes	P280
3 & 5 MHz probes	P350
3 & 8 MHz probes	P380
5 & 8 MHz probes	P580
2 MHz waterproof probe	P2WP

Technical Specifications

Weight:	0.285 kg (10 oz.)
Dimensions:	3.17 x 5.71 x 10.79 cm (1.25 x 2.25 x 4.25 in.)
Power Source:	110V or 220V recharger 2 replaceable AA NIMH batteries
Warranty:	5 years parts 1 year labor
This product is not CE marked	



Nicolet
VASCULAR
a division of VIASYS Healthcare

Medical Device Depot Inc.
www.medicaldevicedepot.com
1-877-646-3300

VIASYS™
HEALTHCARE
Excellence For Life

CE 0344
ISO 9001 CERTIFIED