

#2587/#2588 Power-Lift Designer Blood-Drawing Chairs



ISO 9001-2000 CERTIFIED

Adjustable Height for Patient & Staff Convenience at the Touch of a Button!

Patented blood drawing chair raises and lowers in seconds with portable foot switch. Lower the seat down for safe patient entry and exit; raise the seat up to reduce back strain on the phlebotomist. Electrical components are UL approved for medical use.



#2587
Fiberglass Seat



#2588
Fully Padded Vinyl

U.S. Patent 7,458,638



Proudly Made in the USA!

1-800-237-3377 5516 SW 1st Lane Ocala, FL 34474

www.wincomfg.com

#2587/2588 Power-Lift Designer Blood-Drawing Chairs



ISO 9001-2000 CERTIFIED



W/ Fiberglass seat also available
MODEL # 2587

Electrical components are UL approved for medical use

VINYL SELECTION

(Printed vinyl colors may vary slightly from actual colors)

(A)	02 Royal Blue	03 Taupe	04 Moss Green	06 Hunter Green	07 Grey	08 Black	14 Burgundy	16 Mauve	17 Blueridge			
(B)	32 Grotto	33 Eucalyptus	37 Navy	60 Pyrite	61 Gypsum	62 Ice Blue	63 Cypress					
(C)	13 Granny Smith	25 Spanish Rose	26 Turmeric	27 Steel Grey	28 Carob	29 Curry	30 Doe	31 Tallow	39 Toffee	65 Sage	66 Current	67 Periwinkle


SPECIFICATIONS: #2587

SEAT HEIGHT:	22.5"-28.5" (57.2 cm-72.4 cm)
SEAT WIDTH:	18.5" (47 cm)
OVERALL HEIGHT:	38"-43.5" (96.5 cm-110.5 cm)
ARMREST HEIGHT:	30"-41.5" (76.2 cm-105 cm)
OVERALL WIDTH:	27.5" (69.9 cm) w/o Side Table 38" (96.5 cm) w/ Side Table
WEIGHT:	64-72 LBS. (29-32 kg)
WEIGHT CAPACITY:	300 LBS. (136 kg)
FOOTPRINT:	27.5" X 25.5" (69.9 cm X 64.8 cm)
FOOTPRINT W/ TABLE:	38" X 25.5" (96.5 cm X 64.8 cm)
PENDANT CORD LENGTH:	72" (182.88 cm) FROM REAR OF CHAIR
POWER CORD LENGTH:	120" (304.8 cm) FROM REAR OF CHAIR
*Weight is shown as a max-minimum range and varies depending on chair options ordered	

#2588

SEAT HEIGHT:	25.5"-31" (64.8 cm-78.7 cm)
SEAT WIDTH:	19" (48.3CM)
OVERALL HEIGHT:	47"-52.5" (119.4CM-133.4CM)
ARMREST HEIGHT:	31"-42.5" (78.7CM-108CM)
OVERALL WIDTH:	27.5" (69.9CM) w/o Side Table 38" (96.5CM) w/ Side Table
WEIGHT:	79-87 LBS. (35.8-39.5KG)
WEIGHT CAPACITY:	300 LBS. (136 KG)
FOOTPRINT:	27.5" X 27" (69.9CM X 68.6CM)
FOOTPRINT W/ TABLE:	38" X 27" (96.5CM-68.6CM)
PENDANT CORD LENGTH:	72" (182.88 cm) FROM REAR OF CHAIR
POWER CORD LENGTH:	120" (304.8 cm) FROM REAR OF CHAIR

OPTIONS



TB-133
TB