

PC-Based Cardio Resting ECG

Providing patient state-of-the-art care in resting environments—we care about you

and your patients. Nasiff Associates, Inc. is the first and continues to take the lead in PC-Based Resting ECG systems. Reliable and user-friendly. Our ECG systems are compact and durable. Performing fast testing and managing them easily with you at the point-of-care.

The Cardio Resting ECG offers a practical software and database management system with training, outstanding customer support and legendary quality.




(sample screen)

CardioCard: PC-Based Resting ECG*/PC EKG System**

(for desk-top and notebooks)

- Simultaneous Interpretive 12-Lead Acquisition of your PC ECG/PC EKG
- ECG Interpretation Classes: MIs, Blocks, Enlargements and Axis
- WINDOWS: Easy Point and Click User Interface
- Real-Time 1, 2, 3, 6 or 12- ECG Lead Color Display of ECG/EKG Complex with Real Time ST and HR
- 10-20 ft ECG Patient Cable w/Locking Replaceable Leads
- Full Patient Demographics and Information in an Integrated and Practically Unlimited Cardio Patient Database



- HR w/Global Measurements, Complete Interval Values and ST Information
- Capture Unlimited ECG/EKG Durations of Patient ECG/EKG Data to Hard Drive, Floppy or any Storage Medium
- USB Connection (no batteries or extra power needed)

- User Selectable Lead Formats on Final ECG/EKG Report
- Print Final Report to any Windows Printer (no fading, low cost, readily available)
- ECG to PDF, FAX and E-mail Capability
- Networking Built-in
- ECG/EKG Serial Historical Comparisons



* Resting PC ECG
** Electrocardiograph

- Pacemaker Detection Option
- ECG Measurement Resolutions: 2-4 msec, 1.22uV at 10mm/sec gain
- ECG Frequency Response: 0.05-150Hz std (higher available)
- Heart Rate Range: 30-300bpm
- Lead Leakage: < 10uA
- Input Impedence > 100Mohm
- Gains: 2.5, 5, 10, 10/5, 20 mm/mV CMRR: 120db
- 12Lead Sampling Rate: 250-1000Hz
- A/D Rate: 10kHz
- Time Base: 25 and 50 mm/sec
- Power Consumption: < 200mA at 5V from PC Bus

The Cardio Database Management System is included free with all of our products! This integrates histories for all patients, whether stress, resting, blood pressure, or holter. We are currently developing new tests which we will give you more details on soon.

Medical Device Depot Inc.
www.medicaldevice depot.com
1-877-646-3300

Nasiff Associates, Inc.
leading patient care