



# SELECT SERIES TABLES



**Medical Device Depot**  
Inc.

www.medicaldevicedepot.com  
1-877-646-3300



## TREATMENT TABLE FEATURES

- Easy to clean laminate surfaces
- Vertical spanner adds support to panel legs
- 2" thick foam
- Welt strip conceals upholstery staples
- Scratch and mar resistant surfaces
- 350 lbs. capacity
- Chip resistant Slate Gray, PVC edge trim
- All stress points reinforced by steel
- 15 laminate colors
- 16 vinyl colors
- Optional color graphics (see pages 10 and 11)
- 3-year limited warranty
- Full line of accessories available (see back cover)

## SCALE TABLE FEATURES

- Comfortable work height of 36"
- Top is rimmed on three sides for safety
- Top's front and back rim edges are beveled and the corners rounded for safety
- Easy-clean, replaceable center pad
- Built-in infantometer with folding ends for measuring height
- Easy access to weight, measurement and storage
- Built-in concealed paper dispenser and cutter
- Recessed front toe kick
- Adjustable built-in levelers provide a stable surface
- Easy-clean all laminate base
- Optional storage area locks available
- Ample storage with various door and drawer configurations
- 4 tables to choose from

## SCALE FEATURES


- Dependable, accurate Tanita scale
- Easy-read digital display
- Auto calibration on start-up
- Auto shut-off
- Electric or battery operation (6 AA batteries required)
- Weights in increments of lbs., ozs. or kilos
- Maximum weight 40 lbs.



6 Top Rimmed on 3 Sides



Beveled Sides



Infantometer



Digital Readout

# SELECT SERIES SCALE TABLES

	Length	Height	Width
<b>7820</b>	43"	36"	24"

**Pediatric Scale/Treatment Table**

- Features 2 doors with 1 adjustable shelf



	Length	Height	Width
<b>7830</b>	43"	36"	24"

**Pediatric Scale/Treatment Table**

- Features 2 doors and 1 drawer storage



	Length	Height	Width
<b>7810</b>	43"	36"	24"

**Pediatric Scale/Treatment Table**

- Features 1 drawer with 2 open shelf storage areas



	Length	Height	Width
<b>7840</b>	43"	36"	24"

**Pediatric Scale/Treatment Table**

- Features 2 doors and 2 drawer storage

# SELECT SERIES TREATMENT TABLES



	Length	Height	Width
<b>7133</b>	72"	31"	27"

### Treatment Table

- Features adjustable backrest, storage shelf and 2 Euro-style drawers

	Length	Height	Width
<b>7130</b>	72"	31"	27"

### Treatment Table

- Same as above without adjustable backrest

	Length	Height	Width
<b>7174</b>	72"	31"	27"

### Cabinet Style Laminate Treatment Table

- Features 4 doors, adjustable backrest and 2 storage compartments each with an adjustable shelf

	Length	Height	Width
<b>7178</b>	72"	31"	27"

### Cabinet Style Laminate Treatment Table

- Same as 7174 without backrest



	Length	Height	Width
<b>7114</b>	46"	36"	24"

### Infant Treatment/Changing Table

- Features full shelf, rimmed top on 3 sides, removable pad and 2 round stainless steel top rails

# SELECT SERIES TREATMENT TABLES



	Length	Height	Width
<b>7122</b>	72"	31"	27"

### Straight Line Treatment Table

- Features full length storage shelf and adjustable backrest

	Length	Height	Width
<b>7120</b>	72"	31"	27"

### Straight Line Treatment Table

- Same as above without adjustable backrest

	Length	Height	Width
<b>7144</b>	72"	31"	27"

### Cabinet Style Laminate Treatment Table

- Features an adjustable backrest, 2 storage compartments with sliding doors and an adjustable shelf in each compartment

	Length	Height	Width
<b>7140</b>	72"	31"	27"

### Cabinet Style Treatment Table

- Same as above without adjustable backrest



	Length	Height	Width
<b>7177</b>	72"	31"	27"

### Cabinet Style Treatment Table

- Features 2 sliding doors front and back plus adjustable backrest
- Shipped fully assembled

	Length	Height	Width
<b>7170</b>	72"	31"	27"

### Cabinet Style Treatment Table

- Same as above without adjustable backrest