

STAINLESS STEEL PRODUCTS



	Length	Height	Width
SS-150	14 1/4"	9"	11 1/4"

Stainless Steel Step Stool with Hand Rail

- All welded construction
- Rubber tread top for safety
- 7/8" diameter handrail
- Rubber feet

	Length	Height	Width
SS-140	14 1/4"	9"	11 1/4"

Stainless Steel Step Stool

- All welded construction
- Rubber tread top for safety
- Rubber feet

Medical Device Depot Inc.

www.medicaldevicedepot.com
Toll Free 1-877-646-3300



NOTE: The above stainless steel items are not to be used near an M.R.I.

	Diameter	Total Height
SS-160	14 1/2"	12"

Stainless Steel Kick Bucket and Frame

- Heavy duty stainless steel frame
- 360° bumper on frame
- 2" rubber wheel ball bearing casters
- 13 quart stainless steel bucket

	Diameter
SS-166	14 1/2"

Stainless Steel Kick Bucket Frame Only



	Seat Diameter	Height Range
SS-2142	14"	25"-29 1/2"

Stainless Steel Stool

- All welded frame
- Smooth machine screw height adjustment
- 2" rubber wheel ball bearing casters
- 4" Thick seat for comfort
- Round foot ring
- Easy clean upholstery

	Seat Diameter	Height Range
SS-2142-21	14"	25"-29 1/2"

Stainless Steel Stool with Backrest

- Same as SS-2142 with padded backrest

	Base Diameter	Height Range
SS-31	22"	60"-103"

Stainless Steel IV Pole

- Stainless steel construction
- 2 ram's horn hooks
- 2" rubber wheel ball bearing casters

	Base Diameter	Height Range
SS-21	18" x 22"	31"-50"

Stainless Steel Instrument Stand

- Stainless steel frame
- Removable stainless steel tray
- U-shaped base
- 2" rubber wheel ball bearing casters
- Single knob adjustment