

STANDARD LAB SERIES

FEATURES

- 1 1/4" square, heavy duty, all-welded, tubular steel frame
- Adjustable height and depth flip arm
- Armrest locks in place when in use, then flips out of the way for easy exit
- Neutral Grey powder-coated frame
- 2 knobs lock chrome armrest into position
- Adjustable foot levelers
- Complete series of chairs to choose from
- 300 lbs. capacity



One-piece flip-arm with edge-bumper stripping



2 knobs lock chrome armrest into position

	Seat	Height
6029-P	17" x 16"	20"
	Width	Depth
Overall	43"	26"
	Height Adjustment	
Arms	27"—33"	

Blood Drawing Chair with Flip Arm and Drawer

- Easy-clean, one-piece plastic seat
- Adjustable laminate flip-arm
- Easy-clean one piece drawer
- Optional padded arms



Medical Device Depot Inc.

www.medicaldevicedepot.com
1-877-646-3300

STANDARD LAB SERIES



	Seat	Height
6060-P	17" x 16"	20"
	Width	Depth
Overall	34"	26"
	Height Adjustment	
Arms	27"—33"	

Blood Drawing Chair with Flip Arm

- Easy-clean, one-piece plastic seat
- Adjustable laminate flip-arm
- Optional padded arms

Arm in up position



	Seat	Height
6011-P	17" x 16"	26"
	Width	Depth
Overall	34"	26"
	Height Adjustment	
Arms	34"—40"	

Extra-Tall Blood Drawing Chair with Flip-Arm

- Easy-clean, one-piece plastic seat
- Adjustable laminate flip-arm
- Optional padded arms



	Seat	Height
6022-P	17" x 16"	26"
	Width	Depth
Overall	43"	26"
	Height Adjustment	
Arms	34"—40"	

Extra-Tall Blood Drawing Chair with Flip-Arm and Drawer

- Easy-clean, one-piece plastic seat
- Adjustable laminate flip-arm
- Easy-clean, one-piece drawer
- Optional padded arms



675	Option
Optional Padded Flip-Arm for Standard Lab Series	
<ul style="list-style-type: none"> ▪ 2" padded armrests ▪ Please specify color when ordering 	