

# STANDARD LAB SERIES

## FEATURES

- 1 1/4" square, heavy duty, all-welded, tubular steel frame
- Adjustable height and depth flip arm
- Armrest locks in place when in use, then flips out of the way for easy exit
- Neutral Grey powder-coated frame
- 2 knobs lock chrome armrest into position
- Adjustable foot levelers
- Complete series of chairs to choose from
- 300 lbs. capacity



One-piece flip-arm with edge-bumper stripping



2 knobs lock chrome armrest into position

	Seat	Height
<b>6029-P</b>	17" x 16"	20"
	Width	Depth
Overall	43"	26"
	Height Adjustment	
Arms	27"—33"	

### Blood Drawing Chair with Flip Arm and Drawer

- Easy-clean, one-piece plastic seat
- Adjustable laminate flip-arm
- Easy-clean one piece drawer
- Optional padded arms



**Medical Device Depot** Inc.

www.medicaldevicedepot.com  
Toll Free 1-877-646-3300

# STANDARD LAB SERIES



	Seat	Height
<b>6060-P</b>	17" x 16"	20"
	Width	Depth
Overall	34"	26"
	Height Adjustment	
Arms	27"—33"	

### Blood Drawing Chair with Flip Arm

- Easy-clean, one-piece plastic seat
- Adjustable laminate flip-arm
- Optional padded arms

Arm in up position



	Seat	Height
<b>6011-P</b>	17" x 16"	26"
	Width	Depth
Overall	34"	26"
	Height Adjustment	
Arms	34"—40"	

### Extra-Tall Blood Drawing Chair with Flip-Arm

- Easy-clean, one-piece plastic seat
- Adjustable laminate flip-arm
- Optional padded arms



	Seat	Height
<b>6022-P</b>	17" x 16"	26"
	Width	Depth
Overall	43"	26"
	Height Adjustment	
Arms	34"—40"	

### Extra-Tall Blood Drawing Chair with Flip-Arm and Drawer

- Easy-clean, one-piece plastic seat
- Adjustable laminate flip-arm
- Easy-clean, one-piece drawer
- Optional padded arms



<b>675</b>	<b>Option</b>
<b>Optional Padded Flip-Arm for Standard Lab Series</b>	
<ul style="list-style-type: none"> <li>▪ 2" padded armrests</li> <li>▪ Please specify color when ordering</li> </ul>	