

TruTEMP

NO-TOUCH THERMOMETER

No-Touch – Accurate and Dependable Readings

The TruTEMP no-touch thermometer delivers precise accurate readings, time after time, without the need for disposable probes. The quick and easy-to-use design offers a one-click-to-read temperature measurement in less than a half second – simply point the TruTEMP at the forehead.



Thermal Energy Technology

The TruTEMP No-Touch measures the highest temperature on the forehead. From this value, the TruTEMP estimates the T_{core} or Core Temperature of the body using a proprietary algorithm that compensates for ambient temperature, delivering a precise, accurate reading.

TruTEMP Features:

- Three year warranty
- On/Off audible readout in 3 languages
- Small, convenient, easy-to-use
- Thermal energy technology eliminates patient contact
- Stores previous 32 measurements
- Convenient one-click operation

*Ellie is Ellie Block, SDI's star Customer Service emeritus.

Tru TEMP™

TruTEMP offers an easy-to-read, bright backlit LCD display, clearly visible even in dimly-lit surroundings. Unlike other thermal energy devices that give only skin surface temperatures, TruTEMP calculates core body temperature off the forehead, giving a fast, accurate response.



No Touch
Makes Life Easy



No need to disturb sleeping babies or bed-ridden patients. TruTemp never contacts the patient. Precise, accurate temperature readings are easily measured.



Technical Specifications

Operating Temperature: 50° F - 140° F (10° C - 40° C)

Humidity: ≤ 85%

Power: 3 VDC (2 AA batteries)

Size: 6.1 x 3.5 x 1.6 in (L x W x H)

Weight: 6.2 oz.

Display Resolution: 0.1° F (0.1° C)

Measuring Range:

In body mode: 90° F - 109° F (32° C - 42.9° C)

In Surface Temp mode: 32° F to 244.4° F (0° C - 118° C)

Back-light Display

Precision: From 96.8° F to 102.2° F (36° C to 39° C)
= ± 0.2° C/F

Consumption: ≤ 100 mW

Accuracy: ± 0.54° F

Measuring distance: 2 - 3 in.

Automatic turn-off: 5 sec.

TruTemp Part No. 29-6100

Your dealer is:

Medical Device Depot Inc.

www.medicaldevicedepot.com

1-877-646-3300