

# VALUE SERIES

## FEATURES

- 1" square, all-welded steel frame
- 360° non-locking swing arm
- Height adjustable armrests
- All laminate arms with tee-molded edge bumpers
- Easy-clean, one-piece plastic seat
- Optional padded armrests
- Neutral Grey powder-coated frame
- 2 knobs lock chrome armrest into position
- Adjustable levelers
- 300 lbs. capacity



	Seat	Height
<b>6029-V</b>	17" x 16"	20"
	Width	Depth
Overall	43"	26"
	Height Adjustment	
Arms	27"—33"	

### Blood Drawing Chair with Drawer

- One-piece plastic seat
- Adjustable arms
- Laminate pivot armrest
- Easy-clean, one-piece drawer

	Seat	Height
<b>6060-V</b>	17" x 16"	20"
	Width	Depth
Overall	34"	26"
	Height Adjustment	
Arms	27"—33"	

### Blood Drawing Chair

- One-piece plastic seat
- Adjustable arms
- Laminate pivot armrest





	Seat	Height
<b>6011-V</b>	17" x 16"	26"
	Width	Depth
Overall	34"	26"
	Height Adjustment	
Arms	34"—40"	

**Extra-Tall Blood Drawing Chair**

- Reduces bending for staff
- One-piece plastic seat
- Adjustable arms
- Laminate pivot armrest
- Patient foot rest
- Makes it easy for patient to sit down and exit



	Seat	Height
<b>6022-V</b>	17" x 16"	26"
	Width	Depth
Overall	43"	26"
	Height Adjustment	
Arms	34"—40"	

**Extra-Tall Blood Drawing Chair with Drawer**

- Reduces bending for staff
- One-piece plastic seat
- Adjustable arms
- Laminate pivot armrest
- Patient foot rest
- Easy-clean, one-piece plastic drawer
- Makes it easy for patient to sit down and exit



<b>674</b>	Option
------------	--------

**Optional Padded Armrest**

- Same as standard laminate armrest with the added comfort of upholstery
- Please specify color when ordering