



## > ENDOSCOPY

### Disposables

- > Disposable Sigmoidoscope\_\_\_\_\_ 150
- > Disposable Anoscope\_\_\_\_\_ 150

### Sigmoidoscopes

- > Standard Sigmoidoscope\_\_\_\_\_ 151
- > Fiber Optic Sigmoidoscope\_\_\_\_\_ 151
- > Fiber Optic Sigmoidoscope Set\_\_\_\_\_ 151

### Anoscopes

- > Operating Anoscope\_\_\_\_\_ 152
- > Standard Anoscope\_\_\_\_\_ 152
- > Rotating Anoscope\_\_\_\_\_ 152
- > Fiber Optic Anoscope\_\_\_\_\_ 153

### Endoscopic Lighting

- > Sigmoidoscope/Anoscope  
Lighting Systems\_\_\_\_\_ 154

**Medical Device Depot**  
Inc.

www.medicaldevicedepot.com  
1-877-646-3300

## 150 Disposables

### ENDOSCOPY



Disposable Sigmoidoscope, shown with 3.5 V Halogen Illumination System (Sold Separately)

- 53130** KleenSpec Disposable Sigmoidoscope
- 36019** Fiber-Optic Light Head
- 73500** Halogen 3.5 V Handle Adapter
- 71000-A** 3.5 V Rechargeable Handle

## Disposable Sigmoidoscope

- > Individually wrapped for safety and convenience
- > Highly visible graduations
- > 19 mm x 25 cm

**53130** KleenSpec® Disposable Sigmoidoscope with Obturator (25/box, 4 boxes/case)

### Disposable Sigmoidoscope Accessories

- 36103** Complete Illumination System for Disposable Anoscopes (Includes: 73210, 36019, 73305 and 30200 Insufflation Bulb)
- 73220** 6 V Light Handle and Transformer
- 73210** Light Handle with 3 ft Cord
- 73305** Transformer with 5 ft Cord
- 36019** Fiber Optic Light Head
- 73500** 3.5 V Halogen Handle Adapter
- 71000-A** 3.5 V Direct Plug-in Rechargeable Handle



**36019** Fiber-Optic Light Head

**73210** 6 V Light Handle

**53110** KleenSpec Disposable Anoscope

**73305** 6 V Transformer

## Disposable Anoscope

A high-performance, disposable anoscope for convenient, safe, and efficient examination.

- > 45 degree bevel at distal end to facilitate instrumentation
- > Convenient disposable speculum and obturator ready for instant use
- > 19 mm x 10 cm


**53110** KleenSpec Disposable Anoscope with Obturator (25/box, 4 boxes/case)

### Disposable Anoscope Accessories

- 36103** Complete Illumination System for Disposable Anoscopes (Includes: 73210, 36019, 73305 and 30200 Insufflation Bulb)
- 73220** 6 V Light Handle and Transformer
- 73210** Light Handle with 3 ft Cord
- 73305** Transformer with 5 ft Cord
- 36019** Fiber Optic Light Head



**32820**  
 19 mm x 25 cm  
 Standard Sigmoidoscope,  
 Speculum Only with Obturator



### Fiber Optic Sigmoidoscopes—Standard Sizes

- > Made from autoclavable stainless steel for convenience and durability
- > Hinged window with a neoprene seal avoids fogging during insufflation
- > Graduated tube with removable tip for fiber cleaning

- 32810** 19 mm x 15 cm Standard Sigmoidoscope, Speculum only with Obturator
- 32820** 19 mm x 25 cm Standard Sigmoidoscope, Speculum only with Obturator
- 33830** 19 mm x 30 cm Standard Sigmoidoscope, Speculum only with Obturator
- 32830** 19 mm x 35 cm Standard Sigmoidoscope, Speculum only with Obturator
- 33220** 23 mm x 25 cm Standard Sigmoidoscope, Speculum only with Obturator



**32020**  
 Pediatric/  
 Stricture  
 Speculum Only



### Fiber Optic Sigmoidoscope— Pediatric/Stricture Sizes

- > Made from autoclavable stainless steel for convenience and durability
- > Hinged window with a neoprene seal avoids fogging during insufflation
- > Graduated tube with removable tip for fiber cleaning

- 32010** 11 mm x 10 cm Fiber Optic/Stricture Speculum only with Obturator
- 32020** 11 mm x 25 cm Fiber Optic/Stricture Speculum only with Obturator
- 32410** 15 mm x 15 cm Fiber Optic/Pediatric Speculum only with Obturator
- 32420** 15 mm x 25 cm Pediatric/Stricture Speculum only with Obturator



**35303**  
 Fiber Optic  
 Sigmoidoscope Set

### Fiber Optic Sigmoidoscope Set

**35303 Fiber Optic Sigmoidoscope Set Includes:**

- 32410** 15 mm x 15 cm Pediatric/Stricture Speculum only with Obturator
- 32820** 19 mm x 25 cm Standard Sigmoidoscope, Speculum only with Obturator
- 37023** 23 mm x 7 cm Fiber Optic Anoscope, Speculum only with Obturator
- 73210** Light Handle with 3 ft Cord
- 73305** Transformer with 5 ft Cord
- 30200** Insufflation Bulb, Complete
- 07800-U** 6 V Halogen Lamp
- 30130** 30 cm Suction Tube (standard)
- 05391-U** Fiber Optic Sigmoidoscope Set (35303) Carrying Case

## 152 Standard Anoscopes

ENDOSCOPY



**38614**  
 14 mm Operating  
 Speculum  
 Anoscope

### Operating Anoscope

> Autoclavable for convenience and safety

Anoscope Only 89 mm Speculum Length

<b>38614</b>	14 mm Operating (Cut-out) Speculum Anoscope with Obturator
<b>38619</b>	19 mm Operating (Cut-out) Speculum Anoscope with Obturator
<b>38622</b>	22 mm Operating (Cut-out) Speculum Anoscope with Obturator



**38114**  
 14 mm Standard  
 Speculum  
 Anoscope

### Standard Anoscope

> Autoclavable for convenience and safety

Anoscope Only 89 mm Speculum Length

<b>38108</b>	8 mm Standard Speculum Anoscope with Obturator
<b>38114</b>	14 mm Standard Speculum Anoscope with Obturator
<b>38119</b>	19 mm Standard Speculum Anoscope with Obturator
<b>38122</b>	22 mm Standard Speculum Anoscope with Obturator



**38900**  
 Rotating  
 Anoscope

### Rotating Anoscope

> Autoclavable for convenience and safety

Anoscope Only 68 mm Speculum Length

<b>38900</b>	Rotating Speculum Anoscope with Obturator, 20 mm inside diameter, 22 mm outside diameter
--------------	------------------------------------------------------------------------------------------

8 mm

14 mm

19 mm

22 mm

## Fiber Optic Anoscope

The finest in design and craftsmanship for optimal performance.

- > Autoclavable for convenience and safety
- > 45° bevel tip



**37019**  
 19 mm x 7 cm  
 Fiber Optic  
 Anoscope




**37023**  
 23 mm x 7 cm  
 Fiber Optic  
 Anoscope




**37027**  
 27 mm x 7 cm  
 Fiber Optic  
 Anoscope

### Speculum Only with Obturator


<b>37019</b>	19 mm x 7 cm Fiber Optic Anoscope, Speculum Only with Obturator
<b>37023</b>	23 mm x 7 cm Fiber Optic Anoscope, Speculum Only with Obturator
<b>37027</b>	27 mm x 7 cm Fiber Optic Anoscope, Speculum Only with Obturator



19 mm



23 mm




27 mm

## Sigmoidoscope/Anoscope Lighting Systems

Powerful illumination systems.


> Halogen light provides consistent, long-lasting illumination



## Sigmoidoscope/Anoscope Lighting Systems

Powerful illumination systems.

> Halogen light provides consistent, long-lasting illumination



---

### Sigmoidoscopes/Anoscopes and Lighting Systems Accessories

- > Halogen light provides consistent, long-lasting illumination
  - > Fiber optics provide a cool distal light with no reflections or obstructions (36019, 38700 only)
  - > 73500 adapter fits most Welch Allyn 3.5 V power handles
  - > 73550 adapter fits most Welch Allyn 2.5 V power handles (except PocketScope)
  - > 73220 6 V Light Handle and Transformer provide bright halogen illumination
- 36103** Complete 6 V Illumination System for Disposable Anoscopes (Includes: 36019, 73210, 73305 and 30200)
- 73220** 6 V Light Handle and Transformer (includes 73210 and 73305)
- 73210** 6 V Light Handle only with 3 ft Cord
- 73305** 6 V Transformer only with 5 ft Cord
- 36019** Fiber Optic Light Head for Disposable Sigmoidoscopes and Anoscopes
- 30200** Rectal Insufflation Bulb, Complete
- 30210** Disposable Insufflation Bulb Filters Prevent Contamination of Insufflation Bulb (pack of 10)
- 73500** 3.5 V Halogen Handle Adapter
- 73550** 2.5 V Halogen Handle Adapter
- 38700** Fiber Optic Light Carrier for use with Standard Anoscopes
- 07800-U** 6 V Halogen Replacement Lamp
- 03100-U** 3.5 V Halogen Replacement Lamp
- 03400-U** 2.5 V Halogen Lamp
- 30130** 30 cm Suction Tube (standard)
- 30140** 40 cm Suction Tube