



> THERMOMETRY

Electronic Thermometry

- > SureTemp Plus 692 _____ 100
- > SureTemp Plus 690 _____ 100

Ear Thermometry

- > Braun ThermoScan PRO 4000 _____ 102

100 Electronic Thermometry

THERMOMETRY



01692-200
SureTemp Plus 692
Electronic Thermometer

SureTemp® Plus 692

Easy to use with intuitive operation, reduces the risk of cross-contamination with removable probe well and provides more reliable readings with improved accuracy. Takes approximately 4- to 6-second oral, 10-second pediatric axillary, 15-second adult axillary and 10-second rectal temperatures.

- > Ergonomic design—comfortable to hold, easy to use, and very durable
- > Interchangeable, removable probe wells reduce the risk of cross-contamination
- > Large backlit LCD display with location identification capability
- > Includes a 60-second pulse timer for counting pulse and respiration
- > Mechanical or electronic security options to help prevent theft
- > Takes approximately 5,000 readings on three “AA” batteries

- 01692-200** SureTemp Plus 692 Electronic Thermometer with Wall Mount, Security System with ID Location Field, 4 ft Cord and Oral Probe with Oral Probe Well
- 01692-201** SureTemp Plus 692 Electronic Thermometer with Wall Mount, Security System with ID Location Field, 4 ft Cord and Rectal Probe with Rectal Probe Well
- *01692-300** SureTemp Plus 692 Electronic Thermometer, One Room, with Wall Mount, 9 ft Cord and Oral Probe with Oral Probe Well
- *01692-301** SureTemp Plus 692 Electronic Thermometer, One Room, with Wall Mount, 9 ft Cord and Rectal Probe with Rectal Probe Well
- 01692-700** SureTemp Plus 692 Electronic Thermometer, with Wall Mount, 9 ft Cord and Oral Probe with Oral Probe Well, with Rolling Stand and Holder for Spare Probe & Well

*** Ideal for one-per-room setups with thermometer locked in the wall holder.**



01690-400
SureTemp Plus 690
Electronic Thermometer,
Wall Mount, Oral Probe
with Oral Probe Well

SureTemp® Plus 690

Easy to use with intuitive operation, reduces the risk of cross-contamination with removable probe well and provides more reliable readings with improved accuracy. Takes approximately 4- to 6-second oral temperatures, 10-second pediatric axillary and rectal temperatures, and 15-second adult axillary temperatures.

- > Interchangeable, removable probe wells reduce the risk of cross-contamination
- > Large LCD displays temperatures in Fahrenheit or Celsius
- > Last temperature recall
- > Convenient storage for 25 probe covers
- > Takes approximately 6,000 readings on 3 AA batteries

- 01690-200** SureTemp Plus 690 Electronic Thermometer, 4 ft Cord and Oral Probe with Oral Probe Well
- 01690-201** SureTemp Plus 690 Electronic Thermometer, 4 ft Cord and Rectal Probe with Rectal Probe Well
- *01690-300** SureTemp Plus 690 Electronic Thermometer, One Room, Wall Mount, 9 ft Cord and Oral Probe with Oral Probe Well
- *01690-301** SureTemp Plus 690 Electronic Thermometer, One Room, Wall Mount, 9 ft Cord and Rectal Probe with Rectal Probe Well
- 01690-400** SureTemp Plus 690 Electronic Thermometer, Wall Mount, 4 ft Cord and Oral Probe with Oral Probe Well
- 01690-700** SureTemp Plus 690 Electronic Thermometer, Wall Mount, 9 ft Cord and Oral Probe with Oral Probe Well, with Rolling Stand and Holder for Spare Probe & Well

*** Ideal for one-per-room setups with thermometer locked in the wall holder.**

SureTemp® Plus 692 and 690 Accessories



05031
 Disposable
 Probe Covers

05031-101	Disposable Probe Covers (1,000 covers, packaged 25 per box)	02893-000	Probe and Well Kit, 4 ft Oral
05031-110	Disposable Probe Covers (10,000 covers, packaged 25 per box)	02893-100	Probe and Well Kit, 9 ft Oral
21326-0000	Wall Holder for 690/692	02891-0000	Oral Probe Well, Blue for 690/692
21330-0000	Holder for Spare Probe and Probe Well for 690/692	02891-1000	Rectal Probe Well, Red for 690/692
06138-000	Calibration Key, Assembly for 690/692	08273-000	Thermometer Stand, 690/692 with Holder for Spare Probe and Well
02892-000	Probe and Well Kit, 4 ft Rectal	21333-0000	Bracket, M600 Stand Upgrade, 690/692 with Holder for Spare Probe and Well
02892-100	Probe and Well Kit, 9 ft Rectal	01802-110	9600 Plus Calibration Tester



21330-0000
 Holder for Spare
 Probe and Well



02893-000
 Oral Probe and
 Well Kit, 4 ft



02891-0000
 Oral Probe
 Well



02891-1000
 Rectal Probe
 Well

767 Wall System Add-on

01690-300 SureTemp Plus 690 Electronic Thermometer, One Room, Wall Mount, 9 ft Cord and Oral Probe with Oral Probe Well



76791-2MPX
 767 Integrated Diagnostic
 System

01690-300
 SureTemp Plus 767 Wall
 Thermometer

102 Ear Thermometry

THERMOMETRY



04000-200
 Braun ThermoScan
 PRO 4000 Ear Thermometer

BRAUN *ThermoScan* PRO 4000

Provides the speed and accuracy you need for all types of patients. The exclusive ExacTemp™ technology provides an active user feedback system to ensure proper positioning and improved accuracy.

- > Improved accuracy and repeatability through a preheated probe tip and sensor
- > Probe cover detection system ensures precise infrared heat transfer for accurate temperatures
- > Automatic probe cover eject button for quick and easy disposal of used probe covers
- > Easy-to-read, LCD displays temperatures in Fahrenheit or Celsius
- > Three-year limited warranty

- 04000-200** Braun ThermoScan PRO 4000 Ear Thermometer (with "AA" batteries)
- 04000-300** Braun ThermoScan PRO 4000 Ear Thermometer with Recharging Base Station
- 04000-400** Braun ThermoScan PRO 4000 Ear Thermometer with Roll Stand

Braun ThermoScan® PRO 4000 Accessories

- | | |
|--------------------------------------------------------------------------------------------------------|-----------------------------------------------------------------------|
| 05075-800 Braun ThermoScan PRO 4000 Disposable Probe Covers (800 covers, packaged 20 per box) | 53009-000 Braun ThermoScan PRO 4000 Alkaline Battery |
| 05075-005 Braun ThermoScan PRO 4000 Disposable Probe Covers (5,000 covers, packaged 20 per box) | 53020-0000 Braun ThermoScan PRO 4000 Rechargeable Battery Pack |
| 24001-0000 Braun ThermoScan PRO 4000 Recharging Base Station with Rolling Stand | 80199-000 Replacement Power Cord for PRO 4000 Base Station |
| 24001-1000 Braun ThermoScan PRO 4000 Recharging Base Station and Probe Cover Refill Dispenser | 01802-110 9600 Plus Calibration Tester |

