



Green Series Medical Exam Lights

Introducing the Welch Allyn Green Series™ Medical Exam Lights

Whiter, Brighter and Greener



Medical Device Depot Inc.



www.medicaldevice depot.com

1-877-646-3300



A family of exam lights that can fit any healthcare environment.

Green Series Medical Exam Lights

Welch Allyn Green Series™ Medical Exam Lights are among the first to feature white LEDs that provide a brighter, whiter light than ever before and consume less energy. So now you can see truer tissue color during an exam or procedure while also reducing your current cost of ownership.

GS Medical Exam Lights—Whiter, Brighter and Greener

Whiter and Brighter, *to improve patient exams.*

- 5,500° K white light-emitting diodes (LEDs) as opposed to halogen, for a superior whiter, brighter light
- Truer tissue rendition and uniform spot—no dark or hot spots—for enhanced visualization of the exam area
- Industry-leading lumen performance for a longer life and lower cost of ownership



LED VS. HALOGEN

Greener, *for a lower cost of ownership.*

- No bulb replacement for cost savings and less waste
- Consumes less energy for a reduced cost of ownership
- Intensity controls for the GS 600 and GS 900 and minimal heat dissipation for greater comfort and safety for your patients

Designed to aid in infection control *to help keep your patients safer.*

- Touchless auto on/off controls to reduce your contact with the exam light
- New ergonomic design that is easy to clean
- Use with Welch Allyn disposable sheaths to reduce the risk of infection due to cross contamination



Disposable sheaths available



GS Exam Light IV

- Perfect for the OB/GYN office or for specialties that require very focused light in specific areas
- 3x the typical light output and features a highly flexible light pipe that is twice as long as conventional pipes
- Compact, unobtrusive head that ensures pinpoint accuracy and true tissue color during an exam
- 5-year limited warranty

All lights feature touchless controls

GS 300 General Exam Light

- Ideal for primary care facilities, emergency rooms, outpatient facilities and dermatology
- Perfect combination of intense light output and cool operation, all delivered in a compact design
- Affordable price and rugged performance
- 5-year limited warranty

GS 600 Minor Procedure Light

- Ideal for minor procedures and exams in any healthcare environment
- Features 3 LEDs and a wider head specifically for demanding viewing conditions that require exceptionally good light
- 5-year limited warranty

GS 900 Procedure Light

- Perfect for minor procedures/procedures in the hospital, ambulatory care center, outpatient facility or physician's office
- Features 6 LEDs, outstanding durability, extraordinary maneuverability and unparalleled spot quality
 - Ceiling mount option gives physicians the ability to maximize floor space and direct light wherever it is needed during a procedure
 - Minimal heat dissipation and high color temperature
- 5-year limited warranty

A variety of mounting options available:



Note: Mobile stand will automatically ship with #44400 and #44600